

# 6X Folding Fingerprint Magnifier

A compact, fixed-focus tool designed for quick and portable fingerprint examination. It provides 6× magnification using two optically ground and polished glass lenses for clear viewing, and features a folding metal stand that allows it to be set down for hands-free use. Its durable design folds into a leatherette protective case, making it easy to carry and suitable for both fieldwork and lab settings.



FT332

## SPECIFICATIONS

### Construction

- 6X magnifier lens
- Optically ground

### Materials

- Leatherette protective case (included)
- Folding metal stand
- Polished glass lens

### Performance

- Fixed-focus
- Provide stability when examining evidence
- For use in the field and lab

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI