

# 10X General Purpose Magnifier

A compact forensic inspection tool designed for clear, close-up examination of fingerprints and other small evidence. It features a fixed-focus 10× magnification system using two optically ground and polished glass lenses, which provide sharp, distortion-free viewing. A built-in, battery-powered white light source enhances visibility, making it especially useful in low-light environments such as crime scenes.



FT333

## SPECIFICATIONS

### Construction

- 10X magnifier lens
- Optically ground
- Built-in:
  - Battery-powered push button (included)
  - On / off white light source

### Materials

- Polished glass lens

### Dimensions

- Lens
  - (2) 1" diameter

### Performance

- Fixed-focus
- Provides crystal clear magnification
- Very useful at crime scenes

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI