

10X Illuminated Magnifier

A compact tool designed for close visual inspection of fingerprints and other small evidence. It features a fixed-focus 10× acrylic lens that provides clear magnification for detail work, paired with a built-in light source powered by two “C” batteries to improve visibility in low-light conditions.



FT334

SPECIFICATIONS

Construction

- 10X magnification
- Built-in light source
- Sliding on / off switch for the light
- Case for magnifier (included)
- (2) “C” batteries (included)

Dimensions

- Magnifier: 8” L
- Lens: 1.25” diameter

Materials

- Acrylic lens

Performance

- Fixed-focus
- Used for classifying fingerprints or examining other types of evidence

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI