

10X Tripod Fingerprint Magnifier

A tool designed to help examiners closely inspect fingerprint details. It uses two small, precision-polished biconvex plastic lenses to provide 10× magnification, allowing fine details to be seen clearly. The magnifier is mounted on a tripod base, which stabilizes the lens above the surface being examined and frees the user's hands for tasks like ridge counting or comparison work.



FT335

SPECIFICATIONS

Construction

- (2) 10X magnifier lens
- Tripod mount
- Worm screw height adjustment

Materials

- Biconvex plastic lens

Dimensions

- Lens: (2) 1" diameter

Performance

- Height adjustment allows for precise focusing
- Tripod design allows for ridge counter use

MOPEC
GROUP

Mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI