

Drug Facilitated Sexual Assault Toxicology Collection Kit

A comprehensive Drug Facilitated Sexual Assault Toxicology Collection Kit specimen collection system designed for the proper collection, preservation, and transport of biological samples in cases of suspected drug-facilitated sexual assault.



FT346

SPECIFICATIONS

Construction

- Used when victimizing drugs are suspected
- Collects blood and / or urine
- Kit is comprised of:
 - (1) Pre-sealed Kit Box with chain-of-study and bio-hazard label
 - (1) Kit Instruction Sheet / FDA Insert
 - (1) Victim Information Form Printed on 2-ply NCR paper
 - (1) Inner Box Blood Tube and Urine Bottle Holder
 - (1) 100 ml Urine Specimen Bottle
 - (1) Multi-sample Needle
 - (1) Safety Tube and Needle Holder

- (1) PVP Prep Pad (alcohol-free)
- (2) 10 ml Blood Collection Tubes (gray-stoppered) containing
- (3) Tamper-evident Specimen Security Seals
- (1) Resealable Bag containing a liquid-absorbing sheet
- (1) Tamper-evident Kit Box Shipping Seal

Dimensions

- Pre-sealed kit box
6" L x 4.75" W x 2.5" H
- Multi-sample Needle
20 g x 1.5"

Quantity

- 10 kits

Performance

- Identify drugs in a person's system
- Kit instruction sheets are written for medical personnel and assure proper collection and chain of custody
- After the examination, the kit is handed over to the investigating officer for submission and analysis by the crime laboratory

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI