



Number Wound Mapping Marker, Various Ranges

These Wound Mapping Markers are an adhesive labeling system designed to enhance wound documentation in clinical and forensic photography. Featuring number markers with varying ranges per sheet, it allows users to clearly identify and differentiate multiple wounds within photographs, improving organization and accuracy in medical records.



FT356



FT357



FT358

SPECIFICATIONS

Construction

- Number markers
- 2 arrows per sheet
- Helps in wound documentation by keeping photographed wounds organized in the medical record

Materials

- Skin-friendly adhesive

Quantity

- 100 sheets

Dimensions

- Sheet
6.4 cm x 26.7 cm
- Marker
1.4 cm x 1.4 cm

Performance

- Smaller number ranges give flexibility and help reduce waste
- Helps accelerate assessment times
- Clinical setting
Used to track present admission wounds
- Forensic setting
Keeps evidence organized and tell the story

Variations

Model No.	Description	Number Range
FT356	Number Wound Mapping Marker	1 - 12
FT357	Number Wound Mapping Marker	1 - 14
FT358	Number Wound Mapping Marker	1 - 25

MOPEC
GROUP

mopec
ELEVATING PHOTOGRAPHY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI