


# Adhesive-Backed Fluorescent Photo Scales

An adhesive-backed 2" photo scales that fluoresce under UV light for enhanced visibility when photographing evidence treated with fluorescent powders or dyes. Designed for quick application, they provide an accurate measurement reference for forensic documentation.



FT398

## SPECIFICATIONS

### Construction

- Adhesive-backed
- Available in inches and metric
- Printed in bright orange
- Fluoresces under a UV light source

### Quantity

- 50 per pack

### Dimensions

- 2" (50 mm) L

### Performance

- Easy application and removal to any surface
- Are used when photographing evidence that has been treated with fluorescent powders or dyes

**MOPEC**  
GROUP

**Mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI