


L-Shaped Photomacrographic Scale

Scale is made from laminated, non-glare plastic for accurate evidence photography. Features inch measurements, crosshair circles for parallax correction, and contrasting bars with gray scale for precise measurement and color balance in images.



FT407

SPECIFICATIONS

Construction

- L-shaped scale
- Printed in inches instead of the metric system
- 3 crosshair circles to help in correcting parallax

Materials

- Heavy laminated plastic with a non-glare finish

Options & Accessories

- FS408 fluorescent version

Performance

- Alternating black and white bar measurement of over or underexposed photographs
- 18% gray bars for color corrections

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI