


Shirt Pocket UV Long Wave Lamp

A compact lamp that delivers reliable 365 nm UV illumination in a highly portable design. Perfect for investigators on the move, this lightweight lamp fits easily in a pocket while providing the power needed to reveal fluorescent evidence. With simple battery operation and included accessories, it's a dependable, grab-and-go solution for quick inspections in the field.



FT424

SPECIFICATIONS

Construction

- Portable, pocket-sized lamp
- Made with transistors
- Comes with:
 - Instructions
 - Wrist carrying strap
 - 4 "AA" alkaline batteries

Utilities / Power

- 4 Watt UV lamp
- 365 nm light source

Dimensions

- 6.25" L x 2.125" W x 0.75" H

Performance

- Easy to use for quick inspections
- Portable UV lamp for outside investigations

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI