

# Post-Style Leg Holder for Dissection Table

The Mopec Post-Style Leg Holder provides secure support for lower extremities during dissection procedures. Featuring durable stainless steel construction, it improves specimen stability and integrates with compatible Mopec tables.



## SPECIFICATIONS

### Construction

- Includes 2 leg holders and support posts
- Built to support each leg while working
- Leg holders can be used on any of the 4 corners of the cart

### Materials

- 304 stainless steel construction with a polished No. 4 finish

### Dimensions

- 4.3" L x 5.6" W x 24.9" H

### Compatibility

- HA100
- HA200
- HA210
- HA400
- HB1400
- HB1410
- HB1500
- HB1510
- HB1550
- 600022-05

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI