

# Post-Style Book Holder for Dissection Table

The Mopec Post-Style Book Holder provides a durable, stainless steel elevated platform for reference materials during dissection. It keeps manuals, guides, and notes easily accessible while preserving valuable workspace on compatible tables.



## SPECIFICATIONS

### Construction

- Includes 1 book holder and support post
- Built to support a book upright for users and students to read while working
- Can be used on any of the 4 corners of the cart

### Materials

- 304 stainless steel construction with a polished No. 4 finish

### Dimensions

- 12" L x 8.6" W x 25.5" H
- Useable book-holder area: 12" W x 10" H with 2" lip to hold large books

### Compatibility

- HA100
- HA200
- HA210
- HA400
- HB1400
- HB1410
- HB1500
- HB1510
- HB1550
- 600022-05

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI