

Post-Style Instrument Holder for Dissection Table

The Mopec Post-Style Instrument Holder provides a durable, stainless steel elevated storage solution for surgical and dissection tools. It keeps essential instruments organized, accessible, and off the work surface while integrating seamlessly with compatible tables.



SPECIFICATIONS

Construction

- Includes (1) instrument holder, removable and washable tray, and support post
- Built for placing instruments while working
- Can be used on any of the 4 corners of the cart

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 17" L x 15.6" W x 25.5" H
- Useable tray area: 15.2" L x 10.2" W with a .4" lip to hold tools or catch fluids

Compatibility

- HA100
- HA200
- HA210
- HA400
- HB1400
- HB1410
- HB1500
- HB1510
- HB1550
- 600022-05

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI