

Wall Mounted Autopsy Station, Double Sink

Mopec's Wall-Mounted Autopsy Stations are a durable, stainless steel workstation designed for professional pathology and mortuary environments. It offers a spacious, easy-to-clean work surface with efficient drainage, integrated utilities, and features that support smooth workflow and sanitation. With built-in lighting, storage, waste handling, and compatibility with Mopec carts, this line is reliable for everyday use in demanding settings.



1036-12

SPECIFICATIONS

Construction

- TÜV Certified
- Grid plate work surfaces with 1/2" recessed holes for drainage
- (2) 10' Spray hose assembly that can reach the length of the sink
- Center instrument drawers
- (2) sinks with hot and cold-water capabilities
- (2) cart latch works with most Mopec Carts except those with hydraulic bases
- Wall mounted passive ventilation through Facility-provided HVAC

Utilities / Power

- (4) water-resistant GFCI 120V / 60Hz outlets
- Provide 60-65 psi of water pressure at 4 Gallons per minute for proper operation of hydro aspirator. Refer to drawing for plumbing connection sizes
- Facility-provided Ventilation capable of 400CFM at 0.25 w.g.

Materials

- 14-16 gauge 304 stainless steel with a #4 finish
- Pex piping for plumbing purposes

Dimensions

- 29" L X 132" W X 50" H
- Sink Dimensions: 16" L X 30" W X 6.0" D

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI