

Stainless Steel Embalming Station

Mopec's Autopsy Stations are a durable, stainless steel workstation designed for professional pathology and mortuary environments. It offers a spacious, easy-to-clean work surface with efficient drainage, integrated utilities, and features that support smooth workflow and sanitation. With built-in lighting, storage, and waste handling, this line is reliable for everyday use in demanding settings.



1036-9

SPECIFICATIONS

Construction

- Large sink basin
- Hot and cold water mixing faucet mounted on top deck with high spout for use with embalming machines.
- Deluxe aspirator with reverse flow to relieve suction pressure.
- Cold water spray hose assembly with spray nozzle.
- GFCI duplex receptacle with waterproof cover.
- Access panel

Materials

- 304 stainless steel with No.4 finish

Dimensions

- Overall: 40" D X 32" W X 41" (70" included tall faucet) H
- Sink Basin: 18.6" L X 29.25" W X 8" H

Utilities / Power

- Water: Provide 60–65 PSI water pressure and 4 GPM of flow
- Electrical: Provide (1) circuit in 4" X 4" junction box
- 115V/60Hz / 1Ph / 15A

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI