

Electrostatic Pre-filter for Downflow Workstations

A pre-filter designed for the Reveal Ductless Downflow Workstation that captures airborne particles and aerosols before they reach the main filtration system. It helps extend the life of downstream filters while maintaining efficient airflow and safe workstation operation.



For reference only: Reveal Ductless Downflow Workstation

SPECIFICATIONS

Construction

- (1) Filter provided upon purchase of a Ductless Downflow Workstation
- Electrostatic pre-filter used inside the Ductless Downflow Workstation

Utilities / Power

- 2 options:
120 V / 60 Hz
220 V / 50 Hz

Quantity

- 6 per package

Performance

- Check every 30 days
- First layer of filters that uses static electricity to trap particles in the air before reaching the main HEPA filter
- Helps the workstation last longer and perform more efficiently

Variations

Model No.	Options & Accessories	Size
FS028	FS026 Reveal Ductless Downflow Workstation	24"
FS032	FS030 Reveal Ductless Downflow Workstation	36"

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI