

30" Reveal Fuming Chamber

The Reveal Ductless Automatic Cyanoacrylate Fuming Chamber is a professional forensic tool for safe, consistent fingerprint development. Its durable polypropylene design, LCD controls, and Filtration System ensure precise fuming while protecting users from harmful fumes, with secure, customizable features for reliable evidence processing.



FS034 and FS037

SPECIFICATIONS

Construction

- Includes a filtration system:
 - FS039 Electrostatic Pre-filter
 - FS034: (1) FS035 Carbon Filter
 - FS037: (1) FS038 Carbon Filter
- 2 Hanging rods
- Compact fluorescent lighting
- Centrifugal fan

Materials

- Polypropylene
- Tempered glass

Dimensions

- FS034:
 - External dimensions: 30" W x 48" H x 28" D
 - Internal dimensions: 36" H
- FS037:
 - External dimensions: 30" W x 74.88" H x 28" D
 - Internal dimensions: 57.24" H

Utilities / Power

- 2 options:
 - 120 V / 60 Hz
 - 220 V / 50 Hz

Options & Accessories

- FS034: FS036 Mobile Cart

Performance

- Manage full fuming cycle: humidity, temperature, circulation fan, hot plate, and purge
- Customizable sequences for consistent results
- < 53 noise/dBA/1 meter
- 145 CFM airflow
- 80 FPM face momentum

Variations

Model No.	Description	Type
FS034	30" Reveal Fuming Chamber	Benchtop
FS037	30" Reveal Fuming Chamber	Freestanding

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI