

36" Reveal Ductless Fume Hood

The ductless fume hood provides a safe, efficient workspace with excellent visibility and user-friendly controls. Its advanced filtration system and thoughtful design protect operators while allowing flexible, plug-and-use installation.



FS040

SPECIFICATIONS

Construction

- (2) 15 Watt bulbs
- Centrifugal fan
- (1) Electrostatic Pre-filter
- (1) HEPA Filter
- (1) Carbon Filter
- Manual speed controller
- Hour counter
- Electronic blower switch

Utilities / Power

- 2 options:
120 V / 60 Hz
220 V / 50 Hz

Performance

- < 53 noise/BA/1 meter
- 206 CFM airflow
- 100 FPM face momentum
- Clear back and side panels as well as a double-hinged self-locking clear front sash

Dimensions

- External Dimensions:
36" W x 27" H x 31" D
- Internal Dimensions:
19" H

Materials

- Epoxy coated steel or polypropylene

Options & Accessories

- FS041: Carbon Filter
- FS042: HEPA Filter
- FS043: Electrostatic Pre-filters

Compatibility

- FS044: Mobile Cart

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI