

HEPA Filter for 36" Fume Hood

A high-efficiency particulate air (HEPA) filter designed for the 36" Reveal Ductless Fume Hood to capture fine airborne particles generated during laboratory procedures. As part of the Multiplex™ Filtration System, it helps maintain clean airflow and protects both the operator and environment from particulate contamination.



For reference only: Reveal Ductless Fume Hood

SPECIFICATIONS

Construction

- (1) Filter provided upon purchase of a 36" Reveal Ductless Fume Hood
- A high-efficiency particulate air (HEPA) filter used inside the 36" Reveal Ductless Fume Hood

Utilities / Power

- 2 options:
 - 120 V / 60 Hz
 - 220 V / 50 Hz

Performance

- Check every 3-6 months
- Captures fine particles in the air for user and environmental protection

Options & Accessories

- FS040: 36" Reveal Ductless Fume Hood

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI