

Electrostatic Pre-filters for 36" Fume Hood

A pre-filters designed for the Reveal Ductless Fume Hood to capture larger airborne particles before they reach the main filtration system. They help extend the life of HEPA and carbon filters while maintaining proper airflow and system efficiency. Supplied in a convenient pack of 6 for routine maintenance.



For reference only: Reveal Ductless Fume Hood

SPECIFICATIONS

Construction

- (1) Filter provided upon purchase of a 36" Reveal Ductless Fume Hood
- Electrostatic pre-filter used inside the 36" Reveal Ductless Fume Hood

Quantity

- 6 per package

Utilities / Power

- 2 options:
120 V / 60 Hz
220 V / 50 Hz

Options and Accessories

- FSo40: 36" Reveal Ductless Fume Hood

Performance

- Check every 30 days
- First layer of filters that uses static electricity to trap particles in the air before reaching the main and safety filters
- Can be changed while the unit is operating to prevent accidental exposure to harmful substances

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI