

P-DMAC Paper

The P-DMAC Paper is designed for developing latent fingerprints on porous, heat-sensitive surfaces such as thermal paper and receipts. It reacts with fingerprint residues to produce fluorescent ridge detail when viewed under alternate light sources.



FS470

SPECIFICATIONS

Construction

- Comprised of a set:
 - (1) P-DMAC paper
 - (1) Amber sleeve

Materials

- Paper: P-DMAC

Performance

- 6-month shelf life
To extend longevity of this product, always store in the amber sleeve
- Used for the detection of fingerprints on thermal paper and prints on currency
- Best results usually visualized between 475 nm and 515 nm with an orange barrier filter or glasses

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI