

Dust, Sand, & Dirt Hardener Solution

The Hardener Solution is used to stabilize and harden dust, sand, and dirt in footprint or tire track impressions. It helps prevent deformation caused by the weight of plaster, dental stone, or silicone rubber casting materials.



FS542

SPECIFICATIONS

Construction

- Hardens dust, sand, and dirt
- For footprint or tire track impressions

Performance

- Prevents the impression from deforming
- Shake the container for 10 seconds before pouring into a pump spray bottle

Capacity

- 64oz

Options & Accessories

- FS543: Pump Head Cleaning Solution

MOPEC
GROUP

mopec
ELEVATING METHODOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI