

Basic Plus Turbo Kit

These Turbo Kits are designed specifically for crime scene investigations to enhance the detection of latent evidence while maintaining portability and ease of operation. These kits include everything needed to uncover over 99% of detectable evidence, with multiple configurations available to meet a broad range of investigative requirements and applications.



For reference only: Basic Domestic Violence and Sexual Assault Investigation Kit shown

SPECIFICATIONS

Construction

- Handheld portability
- UltraLite ALS system
- Kit includes:
 - (1) UltraLite ALS Handle
 - (1) BMT Head Assembly
 - (1) Power Adapter Pack
 - (1) AC Power Supply
 - (1) Pair Amber Glasses
 - (1) Vehicle Power Supply
 - (1) 20' Extension Cord
 - (1) UltraLite ALS Case with Foam
 - (1) Battery Charger
 - (3) Battery Packs

Performance

- Battery Charger can charge two Battery Packs at the same time
- Provides everything needed to uncover over 99% of detectable evidence
- Locates volatile substances, such as body fluids, blood splatter, fingerprints, hair, etc.
- Other uses include locating evidence on a person, such as bruising and bite marks

Utilities / Power

- AC Power Supply: 90–240 VAC
- Vehicle Power Supply: 12 VDC

Options & Accessories

- FS855: Basic Turbo Kit
- FS857: Complete Turbo Kit
- FS858: Ultra Turbo Kit

Note:

- Select components of this kit are also sold separately

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI