

Ultra Turbo Kit

These Turbo Kits are designed specifically for crime scene investigations to enhance the detection of latent evidence while maintaining portability and ease of operation. These kits include everything needed to uncover over 99% of detectable evidence, with multiple configurations available to meet a broad range of investigative requirements and applications.



For reference only: Basic Domestic Violence and Sexual Assault Investigation Kit shown

SPECIFICATIONS

Construction

- Handheld portability
- UltraLite ALS system
- Kit includes:
 - (1) UltraLite ALS Handle
 - (1) BMT Head Assembly
 - (1) Power Adapter Pack
 - (1) AC Power Supply
 - (1) Vehicle Power Supply
 - (1) 20' Extension Cord
 - (1) UltraLite ALS Case with Foam Set
 - (1) Battery Charger
 - (3) Battery Packs
 - (5) UV, White Turbo, Yellow, Green, and Red Heads
 - (3) Amber, Red and Yellow Glasses
 - (1) Glasses Bag

Performance

- Battery Charger can charge two Battery Packs at the same time
- Provides everything needed to uncover over 99% of detectable evidence
- Comes in a variety of lights allowing investigators to examine a scene quickly and thoroughly
- Locates volatile substances, such as body fluids, blood splatter, fingerprints, hair, etc.
- Other uses include locating evidence on a person, such as bruising and bite marks

Utilities / Power

- AC Power Supply: 90–240 VAC
- Vehicle Power Supply: 12 VDC

Options & Accessories

- FS855: Basic Turbo Kit
- FS856: Basic Plus Turbo Kit
- FS857: Complete Turbo Kit

Note:

- Select components of this kit are also sold separately

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI