

Glass Slab Fingerprinting Stations

These Glass Slab Fingerprinting Stations are designed for precise control of ink application, helping technicians capture clear prints using a glass inking slab. Available in tabletop and pedestal configurations, they provide consistent performance across different work environments.



FT092



FT095



FT093



FT096

SPECIFICATIONS

Construction

- Includes:
- (1) 0.25" thick glass inking slab
- (1 or 2) Cardholders
- Tabletop:
Divided storage compartment
- Pedestal:
Divided storage compartment
Freestanding base

Performance

- Provides a smooth surface for ink application
- Slab has rounded edges to prevent injury

Materials

- Station: Wood covered in black laminate
- Inking slab: Glass
- Cardholder: Polymer
- Base: Steel

Note:

- Components of these Glass Slab Fingerprinting Stations are also sold separately

Variations

Model No.	Type	Dimensions
FT092	Single Tabletop Glass Slab	28" L x 11" W x 11" H
FT093	Single Pedestal Glass Slab	28" L x 11" W x 37.5" H
FT095	Double Tabletop Glass Slab	33" L x 11" W x 11" H
FT096	Double Pedestal Glass Slab	33" L x 11" W x 37.5" H

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI