

Inking Slab and Backing Plate

These Inking Slabs provide a smooth, reliable surface for consistent ink application, while the Backing Plates are optional and elevate the slab for easier handling and cleaning. The slabs are available in glass or stainless steel and come in two sizes. The plates are available in wood and also come in two sizes.



FT097 and FT098



FT099 and FT100

SPECIFICATIONS

Construction

- Glass slab:
0.25" thick
Provided with non-skid feet
- Stainless steel slab:
Pre-drilled for tabletop mounting
- Wood backing plate:
0.25" thick

Performance

- Slabs:
Provide a controlled surface for applying fingerprint ink
- Backing plate:
Optional
Raises the slab off the station for easy cleaning

Materials

- Slabs: Glass and stainless steel
- Backing Plate: Wood
- Feet: Rubber

Note:

- Slabs and Backing Plates are included in certain Fingerprinting Stations and also sold separately

Variations

Model No.	Description	Dimensions
FT097	Glass Inking Slab	4" X 10"
FT098	Glass Inking Slab	5" X 12"
FT099	Stainless Steel Inking Slab	4" X 10"
FT100	Stainless Steel Inking Slab	5" X 12"
FT101	Wood Riser/Backing Plate	4" X 10"
FT102	Wood Riser/Backing Plate	5" X 12"

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI