

Pre-filters for PCR Prep Stations

A pre-filter designed for the PCR Prep Stations and acts as the first stage of air filtration, capturing larger airborne particles before the air reaches the high-efficiency particulate air (HEPA) filter. This helps extend the life of the HEPA filter while maintaining proper airflow.



SPECIFICATIONS

Construction

- Filter(s) provided upon purchase of a PCR Prep Station
- Pre-filter used inside the PCR Prep Station

Performance

- Captures larger airborne particles before the air reaches the HEPA filter

Variations

Model No.	Description	Size
FT209	Pre-filters for PCR Prep Stations	24"
FT213	Pre-filters for PCR Prep Stations	32"

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI