

48" PCR Prep Station

These Prep Stations are designed to reduce the risk of sample contamination through a clean, controlled environment allowing the operator to amplify RNA and DNA. This is done by delivering continuous positive-pressure HEPA filtered air, the station helps prevent airborne contaminants from entering the work area.



SPECIFICATIONS

Construction

- Electrical cord pass-through port
- Sash position switch and overlap
- UV shelf with integrated pipette holder
- 15 min UV light timer
- Built-in 254 nanometer shortwave lamp
- Safety interlocks
- (2) HEPA Filters
- (2) Pre-filters
- (2) UV Bulbs

Performance

- Provides constant positive pressure HEPA filtered air
- Prevents airborne contaminants from entering the work zone
- Class 100 laminar flow air
For every cubic foot of air, there are fewer than 100 particles that are 0.5 microns or larger in size
- Protects operator from exposure to potential harmful UV radiation
- Supports RNA and DNA amplification

Materials

- Station: Polycarbonate
- Work surface: Polypropylene

Options & Accessories

- Other sizes available for the PCR Prep Station
- FT216: UV Bulbs

Variations

Model No.	Description	Size	Voltage
FT214	PCR Prep Station	48"	110 V
FT215	PCR Prep Station	48"	220 V

MOPEC
GROUP



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI