

Elite Ductless Fume Hood

A ductless fume hood designed to provide chemical fume containment with advanced microprocessor controls, continuous monitoring, and gas phase filtration management.



SPECIFICATIONS

Construction

- Built-in filter and application validation
- Touch screen controller
- Automatic airflow sensor calibration
- Built-in gas detection system
- Spill tray

Dimensions

- FT220 and FT221: 30" W
- FT223 and FT224: 42" W
- FT230 and FT231: 54" W

Materials

- Fiberglass

Performance

- 18-Icon display for simple touch control
- Meets applicable OSHA and ANSI standards
- Conforms to UL, CSA and CE electrical requirements
- LCD readout of airflow in linear feet per minute or meters per second

Utilities / Power

- FT220, FT223, FT230: 110V
- FT221, FT224, FT231: 220V

Options and Accessories

- Safety HEPA filter: FT222, FT225, FT232
- Main HEPA filter: FT226 and FT233
- Safety gas phase filter: FT227 and FT234
- Main gas phase filter: FT228 and FT235
- Pre-Filters: FT229 and FT236

Variations

Model No.	Description	Width	Voltage
FT220	Elite Ductless Fume Hood	30"	110V
FT221	Elite Ductless Fume Hood	30"	220V
FT223	Elite Ductless Fume Hood	42"	110V
FT224	Elite Ductless Fume Hood	42"	220V
FT230	Elite Ductless Fume Hood	54"	110V
FT231	Elite Ductless Fume Hood	54"	220V

MOPEC
GROUP



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI