

# Main HEPA Filter for Elite Ductless Fume Hoods

The Main HEPA Filter for the Elite Ductless Fume Hoods is designed to be the main particulate protection as part of the hood's multi-stage filtration system for applications involving powders and dry particulates. Designed for use alongside the required safety filter and pre-filter.



## SPECIFICATIONS

### Construction

- Option for a type of Main Filter used inside the Elite Ductless Fume Hood
- Used with the Safety HEPA Filter for Elite Ductless Fume Hoods

### Performance

- Use for application involving powders or dry particulates

## Variations

| Model No. | Options & Accessories                          | Size        |
|-----------|------------------------------------------------|-------------|
| FT226     | Main HEPA Filter for Elite Ductless Fume Hoods | 42"         |
| FT233     | Main HEPA Filter for Elite Ductless Fume Hoods | 30" and 54" |

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI