

Pre-filter for Elite Ductless Fume Hoods

The Pre-filter for the Elite Ductless Fume Hoods is designed to be the first stage of filtration in the hood's multi-stage filtration system. Designed for use alongside the required safety filter and main filter. Elite Ductless Fume Hoods can be configured with either HEPA or Gas Phase filtration, and these filter types are interchangeable; however, the Pre-filter should be replaced whenever switching between the two.



SPECIFICATIONS

Construction

- (1) Filter provided upon purchase of an Elite Ductless Fume Hoods
- Used inside the Elite Ductless Fume Hood
- Use with a Safety and Main Filter

Quantity

- 6 per package

Performance

- 2 types of interchangeable filters may be used inside the Elite Ductless Fume Hoods:
 - HEPA Filters
 - Gas Phase Filters
- Pre-filter should be replaced on low airflow and whenever switching between HEPA or Gas Phase Filters

Variations

Model No.	Description	Size
FT229	Pre-filter for Elite Ductless Fume Hoods	42"
FT236	Pre-filter for Elite Ductless Fume Hoods	30" and 54"

MOPEC
GROUP



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI