

Evidence Drying Cabinets (Various Sizes)

These Evidence Drying Cabinets are a specialized laboratory cabinet designed for the secure drying of wet evidence while minimizing the risk of cross-contamination and protecting personnel from exposure to biological materials and odors.



SPECIFICATIONS

Construction

- Locked casters
- Drying rod
- Removable shelves
- One door requires (all sold separately):
 - (2) Pre-Filters
 - (1) HEPA Filter
 - (1) Gas Phase Carbon Filter
- Two door requires (all sold separately):
 - (4) Pre-Filters
 - (2) HEPA Filters
 - (2) Gas Phase Carbon Filters

Utilities / Power

- 110 V

Materials

- Polypropylene
- Door: Steel and epoxy-coated
- Drying rod and shelves: Stainless steel

Options & Accessories

- FT248: Gas Phase Carbon Filter
- FT249: HEPA Filter
- FT250: Pre-filters
- FT251: Waste Pump
- FT252: UV Light

Performance

- 4-stage filtration system: dual pre-filtration, HEPA filter, carbon filter
- Casters for mobility
- 2-speed blower
- Airflow increases when the door(s) are opened
- Provide a secure and dedicated area for drying wet evidence items
- Creates a barrier to keep out the potential threat of cross contamination

Note:

- Pre-filters, HEPA Filters, and Carbon Filters are required and sold separately; Waste Pump and UV Light are optional

Variations

Model No.	Description	Size
FT241	Single Evidence Drying Cabinet	32"
FT242	Single Evidence Drying Cabinet	36"
FT243	Single Evidence Drying Cabinet	48"
FT245	Double Evidence Drying Cabinet	64"
FT246	Single Evidence Drying Cabinet	72"
FT247	Double Evidence Drying Cabinet	72"

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI