

Pre-filters for the Evidence Drying Cabinets

The Pre-filters for the Evidence Drying Cabinets serve as the 1st and 2nd stage of the 4-stage air filtration system. These electrostatically charged filters are designed to capture and remove gross particulate and fibers, helping lower the risk of cross-contamination. These filters should be removed between cases to for examination of any captured trace evidence.



For reference only: shape and size of filter may vary

SPECIFICATIONS

Construction

- Sold separately from purchase of the Evidence Drying Cabinets
- (2) Pre-filters required for Single Evidence Drying Cabinets
- (4) Pre-filters required for Double Evidence Drying Cabinets

Quantity

- 10 per package

Performance

- 1st and 2nd stage of the filtration process for the Evidence Drying Cabinets
- Electrostatically charged
- Remove gross particulate and fibers that could potentially cause cross-contamination
- Collects trace evidence present in the air stream
- Should be removed between uses for laboratory examination

Options & Accessories

- FT241: 32" Single Evidence Drying Cabinet
- FT242: 36" Single Evidence Drying Cabinet
- FT243: 48" Single Evidence Drying Cabinet
- FT245: 64" Double Evidence Drying Cabinet
- FT246: 72" Single Evidence Drying Cabinet
- FT:247: 72" Double Evidence Drying Cabinet

Note:

- Pre-filters are required and sold separately for the Evidence Drying Cabinets

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI