

Pre-filters for the 24" Downflow Latent Print Dusting Workstation

The Pre-filter for the 24" Downflow Latent Print Dusting Workstation serves as the first stage of the filtration system, helping capture larger airborne particles before they reach the main filter. Supports filtration efficiency and prolong the main filters life.



For reference only: 24" Downflow Latent Print Dusting Workstation

SPECIFICATIONS

Construction

- (1) Filter provided upon purchase of a 24" Downflow Latent Print Dusting Workstation
- Pre-filter used inside the 24" Downflow Latent Print Dusting Workstation

Performance

- Helps capture larger airborne particles before reaching the main filter
- Prolongs the main filters life

Quantity

- 10 per quantity

Options & Accessories

- FT253: 24" Downflow Latent Print Dusting Workstation, 110 V
- FT254: 24" Downflow Latent Print Dusting Workstation, 220 V

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI