


Wound Photo Scales Bundle

The Wound Photo Scales are a comprehensive photo documentation system designed to improve the accuracy and consistency of wound photography in clinical and forensic settings. Using a patented, color-coded scale design, the system helps reduce subjectivity in images and standardizes documentation regardless of camera angle or user experience.



FT347

SPECIFICATIONS

Construction

- Color-coded scale
- Bundle of:
 - FT348: Vertical Wound Photo Scales
 - FT349: Wound Photo Scales
 - FT350: Wound Photo Scales
 - FR351: Wound Photo Scales
- FT348:
 - Bars are each 5 cm longer than the one next to it
- FT349, FT350, and FT351:
 - Scale of the ring's inner diameter is provided to assist with ring selection

Materials

- Skin-friendly and latex-free adhesive

Quantity

- FT348: 200 sheets
- FT349: 100 sheets
- FT350: 100 sheets
- FT351: 100 sheets

Performance

- Color-coded to size and adjust to the curvature of the body
- Linear scales support length and width of the wound in the photograph
- A QR code on the sheet that links to an in-service video

Note:

- Each Wound Photo Scale is available for individual purchase

Variations

Model No.	Description	Number of Scales	Area of Scales
FT348	Vertical Wound Photo Scale	1	N/A
FT349	Wound Photo Scales	4	10, 25, 50, 100 sq. cm
FT350	Wound Photo Scales	4	10, 25, 50, 100 sq. cm
FT351	Wound Photo Scales	9	1 - 9 sq. cm

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI