

Vertical Wound Photo Scales

The Wound Photo Scales are a comprehensive photo documentation system designed to improve the accuracy and consistency of wound photography in clinical and forensic settings. Using a patented, color-coded scale design, the system helps reduce subjectivity in images and standardizes documentation regardless of camera angle or user experience.



FT348

SPECIFICATIONS

Construction

- Effective when recording linear wounds
- Includes:
 - (1) In-service flyer
 - (1) Permanent marker
 - (16) Extra note labels
- Color-coded bars are each 5 cm longer than the one next to it
- Has a note label, alignment stick figure, body location tab, and 2 blank tabs
- Blank tabs are provided so you can write down your own body locations
- A QR code on the sheet that links to an in-service video

Materials

- Skin-friendly and latex-free adhesive

Dimensions

- 26.7 cm L x 6 cm W

Quantity

- 200 sheets

Performance

- Disposable
- Eliminates the need for multiple staffing and simplifies the medical recording process
- Color-coded that adjusts to the curvature of the body and helps medical professionals determine the size of the wound

Options & Accessories

- FT347: Wound Photo Scales Bundle

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI