


# Wound Photo Scales (Various Sizes)

The Wound Photo Scales are a comprehensive photo documentation system designed to improve the accuracy and consistency of wound photography in clinical and forensic settings. Using a patented, color-coded scale design, the system helps reduce subjectivity in images and standardizes documentation regardless of camera angle or user experience.



FT349



FT350



FT351

## SPECIFICATIONS

### Construction

- Photo documentation system
- Includes:
  - (1) In-service flyer
  - (1) Permanent marker
  - (16) Extra note labels
- Has multiple linear scales on a sheet
- Color-coded scale of the ring's inner diameter is provided to assist with ring selection
- On the edge of the sheet there is a centimeter ruler to measure the wound
- Has a note label, alignment stick figure, body location tab, and 2 blank tabs

- Blank tabs are provided so you can write down your own body locations
- A QR code on the sheet that links to an in-service video

### Materials

- Skin-friendly and latex-free adhesive

### Dimensions

- Sheet  
18.7 cm x 23.2 cm
- Scales of various sizes

### Quantity

- 100 sheets

### Performance

- Uses a photo scale design to reduce the subjectivity in wound photographs
- System is color-coded to size and adjust to the curvature of the body Standardizes the wound photography regardless of the photographer's skills
- Linear scales support length and width of the wound in the photograph

### Options & Accessories

- FT347: Wound Photo Scales Bundle

## Variations

| Model No. | Description        | Number of Scales | Area of Scales         |
|-----------|--------------------|------------------|------------------------|
| FT349     | Wound Photo Scales | 4                | 10, 25, 50, 100 sq. cm |
| FT350     | Wound Photo Scales | 4                | 10, 25, 50, 100 sq. cm |
| FT351     | Wound Photo Scales | 9                | 1 - 9 sq. cm           |

**MOPEC**  
GROUP

**mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI