

6" Magnetic Scales, Various Colors

These Magnetic Scales are flexible forensic measuring tools designed for accurate evidence photography on metal surfaces. Each scale measures 6 inches in length, comes in various colors, and is constructed from non-glare plastic, ensuring clear visibility and minimal reflection in photographs.



FT387



FT388



FT389

SPECIFICATIONS

Construction

- Magnetic
- Marked in both inches and millimeters
- Crosshair circles at each end
- Available in various colors

Dimensions

- 6" L

Quantity

- 10 per package

Materials

- Non-glare plastic
- Magnetic backing

Performance

- Ideal to use for evidence photography on ferrous surfaces, like file cabinets, motor vehicles, etc.
- Crosshair circles assist in photogrammetry

Variations

Model No.	Description	Printing Color	Base Color
FT387	6" Magnetic Scale	White	Black
FT388	6" Magnetic Scale	Black	White
FT389	6" Magnetic Scale	Black	Yellow

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI