

MegaPictoFramus Fractional Ruler

A large, foldable measuring frame designed for applications requiring oversized reference angles with different colors and sizes available. Made from durable, non-corrosive, and waterproof materials, it features high-visibility Optirectilinear® markings that can be read from distances of up to 40 feet.

FT410 - FT413



SPECIFICATIONS

Construction

- Lightweight
- A line of large folding rules for large frames of reference
- Angles are durable
- Metric measurement reading style
- Printed Optirectilinear® graphics

Materials

- Made of non-corrosive, and waterproof materials

Options & Accessories

- FT409 MegaPictoFramus Ruler Feet
Props the ruler into an upright position

Performance

- Can be easily read from long distances
Up to 40 feet away
- The angles are designed to fold up into a compact to approximately 4'4" x 3.25" x 2"

Variations

Model No.	Description	Dimensions	Color
FT410	MegaPictoFramus Fractional Ruler	4' x 4'	White
FT411	MegaPictoFramus Fractional Ruler	4' x 4'	Yellow
FT412	MegaPictoFramus Fractional Ruler	7' 8" x 7' 8"	White
FT413	MegaPictoFramus Fractional Ruler	7' 8" x 7' 8"	Yellow

MOPEC
GROUP

mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI