

Economy Transport Lift, End/Side Loading

Mopec's lifts are engineered to deliver effortless, smooth, and controlled lifting. Designed for safe maneuverability in mortuaries, laboratories, and medical facilities. These lifts enable operators to transport and elevate loads with precision and confidence.



JD4000



*Tray for reference only

SPECIFICATIONS

Construction

- Concealed Transport Lift
- End and side-loading
- Fits 23" and 27" trays (sold separately)
- For 3 tier racks (when lift and racks are on same floor, may only do two tiers if the rack is in a cooler and the lift is not, check dimensions)
- (4) tray stops
- (4) 5" Full-swivel casters with brakes
- Omni-directional/PVC rollers
- Handheld lift controller
- Includes cover

Materials

- 11GA 304 Stainless steel tubing, various gauges of stainless steel sheet metal
- PVC rollers with HDPE omni-directional rollers
- Vinyl cover

Dimensions

- 81" L x 32.75" W x 14.6 +|- .75" (minimum) – 59.5 +|- .75" (maximum) H
- Cover: 82" L x 33" W x 21" D

Capacity

- 375 lbs (tray weight limit)

Utilities / Power

- Battery-operated electromechanical, scissor lifting system
- 24V power system, (2 Rechargeable deep-cycle 12-volt batteries included) battery cord accepts any standard 120V/1hp/60V outlet

Compatibility

- JZ029: Black Vinyl Cover (included and available individually)
- JC Series 23" and 27" trays
- JZ025: Blue Elastane Cover (sold separately)
- Note: Only compatible with Mopec JC trays, not compatible with Mortech Legacy T3626 trays

Note:

- Tray for reference only. Trays sold separately.

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI