

Cadaver Lift, End Access

Mopec's lifts are engineered to deliver effortless, smooth, and controlled lifting. Designed for safe maneuverability in mortuaries, laboratories, and medical facilities. These lifts enable operators to transport and elevate loads with precision and confidence.



JD5001

SPECIFICATIONS

Construction

- End loading only
- (2) full side rails and two end stop gates
- (4) 5" Full-swivel casters with brakes
- Integrated lift controller with scale display

Capacity

- 500 lbs (tray weight limit)

Materials

- 11GA 304 Stainless steel tubing, various gauges of stainless steel sheet metal
- Galvanized steel rollers
- Stainless steel transfer balls

Dimensions

- 95.2" L x 34.4" W (36.4" including side handle) x 13.4" (minimum) – 76.8" (maximum) H
- 15.5" H minimum for bodies over 300 lbs

Utilities / Power

- Battery-operated electromechanical, scissor lifting system
- 24V power system, (2 Rechargeable deep-cycle 12-volt batteries included) battery cord accepts any standard 120V/1hp/60V outlet

Compatibility

- JC Series 23" and 27" trays
- GA100 32" Legacy tray
- Note: Only compatible with Mopec trays, not compatible with Mortech Legacy trays

Note:

- Trays sold separately.

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI