



Ensemble Pegboard Monitor Mount



SPECIFICATIONS

Construction

- Designed for the Mopec Maestro's Ensemble Pegboard
- Designed to mount a monitor on the pegboard of a grossing station
- Has holes to mount both popular monitor mounting bolt patterns found on the rear of monitors

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 7.0" W x 5.5" H x 1.5" D
- Mounting holes: 100mm x 100 mm and 75mm x 75mm

Weight Capacity

- 25 lbs

Compatibility

- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI