

Raised Dissection Platform

Mopec's Raised Dissection Table provides an elevated, stable work surface designed to support efficient specimen handling and dissection procedures. Its durable construction promotes reliable performance in pathology and laboratory environments.



BC090

SPECIFICATIONS

Construction

- Raised platform with a recessed tabletop
- Creased towards drain for added fluid flow
- Non-marking rubber feet

Capacity

- 25 lbs

Materials

- 304 stainless steel with a polished No. 4 finish
- Rubber feet

Dimensions

- Board Area: 21" L X 14" W X 1.25" H
- Overall: 24" L X 17" W X 15.75" H

Compatibility

- BC002: 16" X 12" White dissection board with feet
- BC005: 16" X 12" Cork dissection board with feet
- BC032: 16" X 12" Light blue dissection board with feet
- BC065: 16" X 12" Dark blue dissection board with feet
- **Note:** Dissection boards sold separately

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI