

Counter Down-Back Draft Workstation

Mopec's BF450 Workstation is a compact, self-contained countertop workstation designed for small laboratories performing dissection, trimming, and staining procedures. Featuring durable stainless-steel construction, an integrated shelf with under-shelf lighting, and Mopec's Dual-Draft Ventilation System, this workstation combines downdraft and backdraft airflow to effectively capture and remove harmful fumes through the facility's ventilation system.



BF450

SPECIFICATIONS

Construction

- Includes (1) BF450 counter down-back draft workstation (facility-vented)
- Includes a shelf with light beneath

Materials

- 304 Stainless steel

Dimensions

- 21" L x 36" W x 19" H

Utilities / Power

- 110/120 V AC/60HZ (USA)
- 10-foot hospital-grade cord
- 18" x 4" rectangular vent duct connection

Performance

- Both downdraft and backdraft functions
- Connects to facility ventilation system
- Removes harmful fumes from the work area through the facility's ventilation system

Options & Accessories

- BF400 (Recirculating version of this workstation)
- BF008 Charcoal filter for BF400
Quantity of 2 required
- BF007 Potassium Permanganate filter for BF400
Quantity of 2 required

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI