


Scissors, Dissection, Sharp/Sharp, Straight

These sharp tipped scissors are designed for precise tissue cuts during post-mortem dissection. The straight, sharp/sharp blades provide clean, controlled cutting for detailed dissection procedures.



SPECIFICATIONS

Construction

- Tips: Sharp/sharp
- Cutting Edge: Straight
- Available in Premium or Standard grades

Materials

- Stainless steel

Dimensions

- AA001 & AA101: 4.5" L
- AA002 & AA102: 5.5" L
- AA003 & AA103: 6.5" L

Variations

Model No.	Description	Grade	Dimensions
AA001	Scissors, Dissection	Standard	4.5" L
AA002	Scissors, Dissection	Standard	5.5" L
AA003	Scissors, Dissection	Standard	6.5" L
AA101	Scissors, Dissection	Premium	4.5" L
AA102	Scissors, Dissection	Premium	5.5" L
AA103	Scissors, Dissection	Premium	6.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI