


Scissors, Dissection, Sharp/Blunt, Straight

These scissors feature one sharp tip, one blunt tip, and a straight cutting edge for controlled dissection. They are ideal for a variety of post-mortem dissection procedures.



SPECIFICATIONS

Construction

- Tips: Sharp/blunt
- Cutting Edge: Straight
- Available in Premium or Standard grades

Materials

- Stainless steel

Dimensions

- AA007 & AA107: 4.5" L
- AA008 & AA108: 5.5" L
- AA009 & AA109: 6.5" L

Variations

Model No.	Description	Grade	Dimensions
AA007	Scissors, Dissection	Standard	4.5" L
AA008	Scissors, Dissection	Standard	5.5" L
AA009	Scissors, Dissection	Standard	6.5" L
AA107	Scissors, Dissection	Premium	4.5" L
AA108	Scissors, Dissection	Premium	5.5" L
AA109	Scissors, Dissection	Premium	6.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI