


Scissors, Dissection, Sharp/Blunt, Curved

These scissors feature one sharp tip, one blunt tip, and a straight cutting edge for controlled dissection. They are ideal for a variety of post-mortem dissection procedures.



SPECIFICATIONS

Construction

- Tips: Sharp/blunt
- Cutting Edge: Curved
- Available in Premium or Standard grades

Materials

- Stainless steel

Dimensions

- AA010 & AA110: 4.5" L
- AA011 & AA111: 5.5" L
- AA012 & AA112: 6.5" L

Variations

Model No.	Description	Grade	Dimensions
AA010	Scissors, Dissection	Standard	4.5" L
AA011	Scissors, Dissection	Standard	5.5" L
AA012	Scissors, Dissection	Standard	6.5" L
AA110	Scissors, Dissection	Premium	4.5" L
AA111	Scissors, Dissection	Premium	5.5" L
AA112	Scissors, Dissection	Premium	6.5" L

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI