

Scissors, Dissection, Blunt/Blunt, Curved

These blunt, curved tip scissors are designed for safe and controlled cutting during post-mortem dissection. The curved blades provide improved access and maneuverability when working with tissue.



SPECIFICATIONS

Construction

- Tips: Blunt/blunt
- Cutting Edge: Curved
- Available in Premium or Standard grades

Materials

- Stainless steel

Dimensions

- AA016 & AA116: 4.5" L
- AA017 & AA117: 5.5" L
- AA018 & AA118: 6.5" L

Variations

Model No.	Description	Grade	Dimensions
AA016	Scissors, Dissection	Standard	4.5" L
AA017	Scissors, Dissection	Standard	5.5" L
AA018	Scissors, Dissection	Standard	6.5" L
AA116	Scissors, Dissection	Premium	4.5" L
AA117	Scissors, Dissection	Premium	5.5" L
AA118	Scissors, Dissection	Premium	6.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI