

Scissors, Dubois Decapitation

These scissors feature slim, blunt, curved tips and thin, delicate shafts for precise post-mortem dissection with minimal damage to surrounding tissues. Available in a variety of lengths and grades to accommodate different dissection needs.



SPECIFICATIONS

Construction

- Tips: Blunt
- Available in Premium or Standard grades

Materials

- Stainless steel

Dimensions

- AA027 & AA127: 10.5" L

Variations

Model No.	Description	Grade	Dimensions
AA027	Scissors, Dubois Decapitation	Standard	10.5" L
AA127	Scissors, Dubois Decapitation	Premium	10.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI