

Scissors, Iris, Curved

These scissors are ideal for small, precise cuts during post-mortem dissection. Their design provides controlled handling for detailed dissection procedures.



SPECIFICATIONS

Construction

- Tips: Sharp points
- Cutting Edge: Curved
- Available in Premium or Standard grades

Materials

- Stainless steel

Dimensions

- AA029 & AA129: 4.5" L

Variations

Model No.	Description	Grade	Dimensions
AA029	Scissors, Iris	Standard	4.5" L
AA129	Scissors, Iris	Premium	4.5" L

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI