

Scissors, Doyen Abdominal

These scissors are ideal for small, precise cuts during post-mortem dissection. Their design provides controlled handling for detailed dissection procedures.



SPECIFICATIONS

Construction

- Tips: Blunt points
- Available in Premium or Standard grades
- Available in Curved or Straight cutting edge

Materials

- Stainless steel

Dimensions

- AA035 & AA036 & AA135 & AA136: 7" L

Variations

Model No.	Description	Grade	Edge
AA035	Scissors, Doyen Abdominal	Standard	Straight
AA036	Scissors, Doyen Abdominal	Standard	Curved
AA135	Scissors, Doyen Abdominal	Premium	Straight
AA136	Scissors, Doyen Abdominal	Premium	Curved

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI